



**LEPAGE**

WWW.LEPAGE.EXPERT

Innovative  
medical solutions



mesotech



# Leadergun

## Fully automatic meso-gun

The only one with anti-leakage system and several doses for each different treatment [CE0476]. It was designed to overpass the limits of all other mesoguns:

### Injection

a single injection (spot)  
burst injection (very fast, 3 injections per second) continuous injection (no stop)

### Dosage

pico dose (0.007 -> 0.05) micro doses (0.013 -> 0.10) macro doses (0.025 -> 0.15)  
Penetration / depth: from 1mm to 13mm Syringe seat - (mL): 1; 2.5; 5; 10 Needles: 4mm; 6mm; 13mm.

### SMT – Step Motor Technology

Extreme precision during the injection and the output of the product and during the extraction of the needle from the skin. Thanks to these features SMT is also used on medical devices in Oncology and Intensive Therapy.

### Automatic functions

Automatic turn off of device after few minutes without manual use.  
Automatic adjustment of penetration depth.  
Automatic detection of the empty syringe and recalling.

### Programming

Syringe size  
Needle size  
Penetration depth  
Injected dose (Pico, Micro, Macro) Injection mode  
Injection frequency  
Permanence of the needle into the skin

Leadergun has logical memory to save custom presets. The integrated software is designed to avoid any risk of error by the operator.

LCD Colour Screen 1,5”  
Brightness and Contrast’ s Adjustment

### Syringe Cover

To avoid eventual touch of the syringe during the treatment.

The cover and the backside of Leadergun have a magnet on its top, which allow the cover’s blockage if it’s turned upside down during treatment.

(According to MD Certificates’ Safety Regulations.)



### Five preset programs

Glow (anti-aging)  
Lipo (cellulite)  
Stretch (stretch-marks)  
PRP  
Botulinum

### Handle Design

Leadergun’s handle has been designed to fit all ergonomics standars. The special TPE material which wraps the handle, allows a perfect grip and a better handling during the treatment.

### Ergonomics

All the components in Leadergun are equally balanced and respond to Ergonomics MSI Requirements

### Battery

Lithium Ion 2200 mAh 7,4 V Autonomy : 36 hours  
Charge time : 2 hours Max  
Battery level of charge on LCD Screen

### Power supply

Universal AC input/ Full range, INPUT: 100 – 240 VAC - 50/60 Mz - 0,8 A



Available in black or white version  
Smart color LCD display